Application

The self-adhesive 3-line labels are available without or with process indicator for steam, ethylene oxide, formaldehyde and hydrogen peroxide sterilization processes. They can be printed with the **gke** labelling devices and user-defined number of labels with the identical content can be generated.

Use of labels in double door sterilizers

For double door sterilizers labels without process indicators may be used. Process indicators are not absolutely necessary because the door into the clean room area only opens, if all sterilization process parameters have been achieved. In this case the labels may be adhered in the clean room area after the sterilization processes have been completed successfully.

Use of labels in single door sterilizers

To differentiate between processed and unprocessed packages in the same room, labels with process indicator should be used.

GMP rules and the Medical Devices Directive (MDD) require validated sterilization processes. Besides validation, routine monitoring of all sterilization processes is mandatory to assure sterility of all packs. Routine monitoring may be achieved by using a **gke Steri-Record®** batch or process monitoring system. Chemical packaging indicators can be used in steam sterilization processes with insufficient air removal characteristics (sterilizers with gravity displacement cycles). The documentation of routine monitoring in the CSSD and/or in the patient documentation is required.

By using the hand labelling device, the 3-line label can be printed on the first line with the initials or number of the responsible person, the sterilizer and batch number as well as the package content number if all packages in the load are of the same content. The 2nd and 3rd line documents the production and expiry date of the goods to be sterilized. One label with the same imprint should be adhered on each pack and onto the documentation sheet. Most sterilizers are equipped with an automatic counter for batch numbers which are generated with every sterilization cycle. The batch number can also be generated manually, e. g. serial number and year. Before sterilized goods are used in the operation theatre the double self-adhesive label will be removed from the pack and adhered onto the patient documentation sheet. This simple method enables a patient related documentation of sterilized goods without any additional manual paperwork.

In case of a secondary infection the labels in the patient documentation can directly refer to the documentation of routine monitoring in the CSSD showing successful sterilization parameters and supply documented evidence.

Product Description

The double self-adhesive 3-line labels (2,5 x 3,1 cm) contain three paper layers and two self-adhesive layers in between. The upper layer of hydrogen peroxide labels is made of plastic foil. One roll with 1" core contains 750 labels. Packages of 2, 4 and 12 rolls are available.

The single self-adhesive 3-line labels (2,5 x 3,1 cm) contain two paper layers and one self-adhesive layer in between. One roll with 1" core contains 1.000 labels. Packages of 12 rolls are available.

The documentation labels contain different process indicators for sterilization processes, see table. Labels that are used after the sterilization process do not require an indicator. The label can be printed by using a special **gke** labelling device.

Performance Characteristics

All **gke** labels with process indicators comply with the standard EN ISO 11140-1 type 1 and provide the information that a pack has passed a sterilization process. The labels do not give any information about the sterility of the goods being sterilized. The glue is specially developed for the adhesive strength of textile, paper and foil packaging material in sterilization processes.

Handling information

- The **gke** labeling devices have to be used to print labels. The labeling devices only differ by printing numeric or alpha-numeric letters in the first line.
- About 10 labels of the roll have to be taken off before inserting the roll into the labeling device. Therefore, about 15 additional labels are added on top of the 750 labels on each roll free of charge.
- Insert the labels into the **gke Steri-Record®** labelling device according to the directions for use which come with the labeller.
- The labeller print head has 3 printing lines. The first line can be used for the following information:
The first three letters can be used for the initials of the person responsible and have to be adjusted with the conversion table shown on the labelling device to select the correct letters. The letters shown will not be printed.
The remaining 9 digits can be used to print the number of the sterilizer and the batch number. All numbers shown will be printed.
The second line contains the sterilization date.
The third line contains the expiry date.
- The indicator changes its colour depending on the type of indicator and sterilization process used. The indicator colour and reference colour are printed on each label.

Storage and Disposal

- For longer periods store all indicators in original package.
- Store indicators before and after sterilization in a dry room at temperatures between 5 - 30°C or 41 - 86°F and a humidity of 5 - 80%.
- The vapour of chemicals especially hydrogen peroxide may destroy the indicator. Therefore, do not store them together with other chemicals.
- The indicators should not be used after expiry date and can be disposed with normal garbage.

For further technical details please contact your local dealer or the **gke** application laboratory. We will assist you with any technical questions. Visit our website www.gke.eu for more information.

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760-052 ED V04 01/2019

Directions for Use

for all **gke Steri-Record®** 3-line labels for batch and patient related documentation

Hand Labeler:

Art.-No.*	Product code	Application
240-820	D-G-NL	First printing line contains 3 alpha-numeric and 9 numeric characters
240-830	D-G-AL	First printing line contains 12 alpha-numeric characters

3-line single self-adhesive labels:

Each roll contains 1.000 labels which are available in the following colours and quantities:

Art.-No.*	Product code	Background colour	Application	Sterilization process	Content (number of rolls)
230-864	D-L-SA-Y	Yellow	<u>Without</u> process indicator	Steam, Ethylene Oxide, Formaldehyde, Dry Heat	12
230-872	C-S-L-1-SA-B	Blue	<u>With</u> type 1 process indicator according EN ISO 11140-1	Steam	

3-line double-self-adhesive labels:

Each roll contains 750 labels with one ink-roll. Labels are available in the following colours and quantities:

Art.-No.*	Product code	Background colour	Application/Standard	Sterilization Process	Content (1 ink-roll and number of rolls)
240-861	D-L-DA-G	Green	<u>Without</u> process indicator	Steam, Ethylene Oxide, Formaldehyde, Dry Heat	12
240-862	D-L-DA-B	Blue			
240-863	D-L-DA-R	Red			
240-864	D-L-DA-Y	Yellow			
240-853	D-L-DA-R	Red			4
240-871	C-S-L-1-DA-G	Green	<u>With</u> type 1 process indicator according EN ISO 11140-1	Steam	12
240-872	C-S-L-1-DA-B	Blue			
240-873	C-S-L-1-DA-R	Red			
240-874	C-S-L-1-DA-Y	Yellow			
240-883	C-S-L-1-DA-R	Red			
242-875	C-E-L-1-DA-P	Purple		Ethylene Oxide	12
242-885					2
243-874	C-F-L-1-DA-Y	Yellow		Formaldehyde	12
243-884					2
244-873	C-V-L-1-DA-R	Red		Hydrogen Peroxide	12
244-883					2

3-line double self-adhesive labels:

Each roll contains 750 labels which are available in the following colours and quantities:

Art.-No.*	Product code	Background colour	Application	Sterilization Process	Content (number of rolls)
000-863	D-L-DA-R-ECO	Red	<u>Without</u> process indicator	Steam, Ethylene Oxide, Formaldehyde, Dry Heat	12
000-873	C-S-L-1-DA-R-ECO		<u>With</u> process indicator	Steam	
000-875	C-E-L-1-DA-P-ECO	Purple		Ethylene Oxide	

(*) To the article numbers a 3-digit letter code is added. The additional letter code refers to the language and/or customized version. It is only added on the outside label; the inside of the pack is identical to the article numbers and the above table.