

Documentation Sheet

Test results with *gke Clean Record®* Cleaning Process Indicators (CPI) in bedpan washers



Hospital:	City:
Department:	Date:
Manufacturer:	Type:
Year of manufacture:	Machine no.:
Indicator Lot no.:	Inspector:

The below documented test result is:

<input type="checkbox"/> a reference result for the above machine. Later on, tests are comparable in the identically constructed machines
<input type="checkbox"/> daily <input type="checkbox"/> weekly routine testing Evaluation: <input type="checkbox"/> Test result = reference → OK <input type="checkbox"/> Test result ≠ reference → not OK

Program:	Lot-No.:
----------	----------

Bedpan					
bottom, inside	on the edge, right	on the edge, left	under the edge, right	under the edge, left	bottom, outside
<i>Please adhere indicator here</i>					
handle	<i>additional test position</i>	<i>additional test position</i>	<i>additional test position</i>	<i>additional test position</i>	<i>additional test position</i>

Urine bottle					
lateral, outside	bottom, outside	bottle neck, inside	lateral, inside	bottom inside	<i>additional test position</i>
<i>Please adhere indicator here</i>					