

Practice/Laboratory/Hospital etc



Documentation sheet for GKE Steri-Record® Stearo-Ampoules

Steam sterilization process	
Batch no.	
Expiry Date	
Ampoule no. for vitality test	
Result vitality test	+ / -
Incubation temperature	55 - 60°C
min./max. incubation time	48 hours/ 7 days

No. (*)	Job	Sterilizer/sterilization process			Incubation time				Result (**) + / -	Responsible Person	
		Date	No.	Batch no.	Start		End			Initials	Signature
					Date	Time	Date	Time			
	Vitality test										

766-004 EN V03 03/2023

(*) If more than one ampoule has to be sterilized in one batch please number them by using a waterproof pen.
 (**) Explanation: (+) Growth with colour change = NON-STERILE
 (--) No growth, no colour change = STERILE