
	Technical Information	730-121-EN		V04
	Generating lot numbers and labeling on GKE packages	Created	11.09.2014	KH
		Changed	31.08.2021	KP
		Checked	31.08.2021	UK
		Released	31.08.2021	UK
File no.: 0.0				

Many GKE products consist of more than one component which are manufactured in separated steps and have separate lot numbers. Inside the package are several components which may have individual lot number, e.g. sealing kit, directions for use, package label, documentation sheet etc.

GKE indicates a lot number on all packages, to be able to trace all components. However, this lot number (LOT) on the outside packaging differs from the lot number printed on the main component (e.g. on the backside). The lot number of the main component is added as "Main Component Lot = MCL" on the outside label. The main component is, for example, chemical or biological indicators.

Example of the content of the package label:



Lot number which is printed on back side of the indicator card or on the label of a biological indicator

If one product contains more than one main component (e.g. start-up kit with PCD and indicator strips) the lot number of the indicators will be stated on the package as additional MCL. This applies for introduction packs consisting of PCD and indicators.

Testsets with different PCDs do not contain a main component lot on the package.